


# 2012 VERSANDTERMINE | SHIPPING SCHEDULE



|           |  |  |  |
|-----------|--|--|--|
| <b>1</b>  |  |  |  |
| <b>1A</b> | Capture-R® Ready-Screen® (4), incl. Cw, Kp(a) Plate  | 0066802<br>0066812   | 24 tests<br>120 tests  |
| <b>1B</b> | Capture-R® Ready-Screen® (I and II) incl. C* Plate   | 0066202<br>0066212   | 48 tests<br>240 tests  |
| <b>1C</b> | Capture-R® Ready-Screen® (Pooled Cells), incl. C* Plate  | 0066203<br>0066213   | 96 tests<br>480 tests  |
| <b>1D</b> | Capture-R® Ready-Screen® (3), incl. C*, Kp(a) Plate<br>Capture-R® Ready-ID® Plate<br>Capture-R® Ready-ID® Extend I Plate<br>Capture-R® Ready-ID® Extend II Plate | 0066804<br>0066814<br>0066204<br>0066214<br>0006454<br>0006455<br>0006456<br>0006457 | 24 tests<br>120 tests<br>6 tests<br>30 tests<br>6 tests<br>30 tests<br>6 tests<br>30 tests |
| <b>1E</b> | Capture-P® Plate<br>Capture-P® Ready-Screen® Plate   | 0066246<br>0066251<br>0066247  | 96 tests<br>480 tests<br>6 tests   |

|          |                                     |         |              |
|----------|-------------------------------------|---------|--------------|
| <b>2</b> |                                     |         |              |
|          | Capture-R® Ready Indicator Cells    | 0006427 | 11.5 ml      |
|          | Referencells® A(1), A(2), B and O   | 0002338 | 10 x 11.5 ml |
|          | Referencells® A(1), B               | 0002345 | 4 x 10 ml    |
|          | Referencells® A(2)                  | 0002342 | 2 x 10 ml    |
|          | Panoscreen® I and II                | 0002380 | 10 ml        |
|          | Panoscreen® I, II and III, incl. C* | 0002377 | 2 x 10 ml    |
|          | Panoscreen® I, II and III           | 0002381 | 3 x 10 ml    |
|          | Hemantigen® (Pooled Cells)          | 0002223 | 10 ml        |
|          | Panocell®-16 Plus                   | 0002332 | 3 x 10 ml    |
|          | Checkcell®                          | 0002224 | 16 x 3 ml    |
|          | Checkcell® Weak                     | 0002225 | 10 ml        |
|          | Complement Control Cells            | 0002227 | 3 x 10 ml    |
|          |                                     | 0007930 | 3 ml         |

|          |                                |         |            |
|----------|--------------------------------|---------|------------|
| <b>3</b> |                                |         |            |
|          | Capture-P® Indicator Red Cells | 0066240 | 11.5 ml    |
|          | Panocell®-10                   | 0003032 | 11 x 3 ml  |
|          | Panocell®-10, Ficin-Treated    | 0002385 | + 1 x 2 ml |
|          | WB corQC™                      | 0066090 | 23 x 3 ml  |
|          |                                |         | + 1 x 2 ml |
|          |                                |         | 2 sets of  |
|          |                                |         | 4 x 4.5 ml |

|          |              |         |           |
|----------|--------------|---------|-----------|
| <b>4</b> |              |         |           |
|          | Panocell®-20 | 0005020 | 20 x 3 ml |

|          |                                    |         |             |
|----------|------------------------------------|---------|-------------|
| <b>6</b> |                                    |         |             |
|          | corQC™                             | 0002410 | 2 x 11.5 ml |
|          | corQC™ EXTEND Cell 1, 2 and 3      | 0002417 | 3 x 5 ml    |
|          | corQC™ EXTEND Complete             | 0066297 | 1 x 10 ml   |
|          |                                    |         | + 4 x 5 ml  |
|          | corQC™ EXTEND Standard             | 0066296 | 10 ml       |
|          | DAT Positive Control Cells Galileo | 0066122 | 10 ml       |
|          |                                    | 0066125 | 4 x 10 ml   |

|          |                            |         |      |
|----------|----------------------------|---------|------|
| <b>7</b> |                            |         |      |
|          | Weak D Cells               | 0002995 | 5 ml |
|          | Tech-Chek™                 | 0002999 | kit  |
|          | Fetal Bleed Screening Test | 0007892 | kit  |

|          |                               |         |       |
|----------|-------------------------------|---------|-------|
| <b>8</b> |                               |         |       |
|          | GammaZyme-F™ (Ficin solution) | 0007056 | 10 ml |

|          |       |         |     |
|----------|-------|---------|-----|
| <b>9</b> |       |         |     |
|          | RISE™ | 0007999 | kit |

Auftragsannahme / Customer Service  
 Tel.: +49 (0)6074 8420-20  
 order.germany@immucor.com